



TM158 rabbit pAb

Cat No.:ES19868

For research use only

Overview

Product Name	TM158 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TM158 AA range: 97-147
Specificity	This antibody detects endogenous levels of human TM158
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Transmembrane protein 158 (40 kDa BINP-binding protein) (p40BBP) (Ras-induced senescence protein 1)
Gene Name	TMEM158 HBBP RIS1
Cellular localization	Membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	25907
Human Swiss-Prot Number	Q8WZ71
Alternative Names	
Background	

