



HEAT2 rabbit pAb

Cat No.:ES19881

For research use only

Overview

Product Name	HEAT2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human HEAT2 AA range: 242-292
Specificity	This antibody detects endogenous levels of human HEAT2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	HEAT repeat-containing protein 2
Gene Name	HEATR2
Cellular localization	Cytoplasm . Dynein axonemal particle . Observed only in the cytoplasm of ciliated cells and absent from cilia. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	54919
Human Swiss-Prot Number	Q86Y56
Alternative Names	
Background	

