

SDIM1 rabbit pAb

Cat No.:ES19887

For research use only

Overview

Product Name	SDIM1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SDIM1 AA range: 83-133
Specificity	This antibody detects endogenous levels of human SDIM1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Stress-responsive DNAJB4-interacting membrane protein 1
Gene Name	SDIM1
Cellular localization	Membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q6ZPB5
Alternative Names	
Background	

