

# ECM29 rabbit pAb

Cat No.:ES19892

For research use only

## Overview

Product Name	ECM29 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ECM29 AA range: 1599-1649
Specificity	This antibody detects endogenous levels of human ECM29
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Proteasome-associated protein ECM29 homolog (Ecm29)
Gene Name	ECM29 KIAA0368
Cellular localization	Endoplasmic reticulum . Endoplasmic reticulum-Golgi intermediate compartment . Endosome . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Nucleus . Endosome, multivesicular body . Cytoplasmic vesicle .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q5VYK3
Alternative Names	
Background	

