



PSMG1 rabbit pAb

Cat No.:ES19894

For research use only

Overview

Product Name	PSMG1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PSMG1 AA range: 144-194
Specificity	This antibody detects endogenous levels of human PSMG1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Proteasome assembly chaperone 1 (PAC-1) (Chromosome 21 leucine-rich protein) (C21-LRP) (Down syndrome critical region protein 2)
Gene Name	PSMG1 C21LRP DSCR2 PAC1
Cellular localization	Cytoplasm. Endoplasmic reticulum.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	8624
Human Swiss-Prot Number	O95456
Alternative Names	
Background	

