



CAR17 rabbit pAb

Cat No.:ES19899

For research use only

Overview

Product Name	CAR17 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CAR17 AA range: 22-72
Specificity	This antibody detects endogenous levels of human CAR17
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Caspase recruitment domain-containing protein 17 (Caspase-1 inhibitor INCA) (Inhibitory caspase recruitment domain protein)
Gene Name	CARD17 INCA
Cellular localization	Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	440068
Human Swiss-Prot Number	Q5XLA6
Alternative Names	
Background	

