

# SEBOX rabbit pAb

Cat No.:ES19900

For research use only

## Overview

Product Name	SEBOX rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SEBOX AA range: 109-159
Specificity	This antibody detects endogenous levels of human SEBOX
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Homeobox protein SEBOX (Homeobox OG-9) (Skin-, embryo-, brain- and oocyte-specific homeobox)
Gene Name	SEBOX OG9X
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	645832
Human Swiss-Prot Number	Q9HB31
Alternative Names	
Background	

