

# CR032 rabbit pAb

Cat No.:ES19905

For research use only

## Overview

Product Name	CR032 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CR032 AA range: 14-64
Specificity	This antibody detects endogenous levels of human CR032
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	UPF0729 protein C18orf32 (Putative NF-kappa-B-activating protein 200)
Gene Name	C18orf32
Cellular localization	Endoplasmic reticulum . Lipid droplet .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	497661
Human Swiss-Prot Number	Q8TCD1
Alternative Names	
Background	

