

EIF3M rabbit pAb

Cat No.:ES19908

For research use only

Overview

Product Name	EIF3M rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human EIF3M AA range: 107-157
Specificity	This antibody detects endogenous levels of human EIF3M
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Eukaryotic translation initiation factor 3 subunit M (eIF3m) (Fetal lung protein B5) (hFL-B5) (PCI domain-containing protein 1)
Gene Name	EIF3M HFLB5 PCID1 GA17 PNAS-125
Cellular localization	Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	10480
Human Swiss-Prot Number	Q7L2H7
Alternative Names	
Background	

