

PAQRB rabbit pAb

Cat No.:ES19909

For research use only

Overview

Product Name	PAQRB rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PAQRB AA range: 38-88
Specificity	This antibody detects endogenous levels of human PAQRB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Monocyte to macrophage differentiation protein (Progesterone and adipoQ receptor family member 11) (Progesterone and adipoQ receptor family member XI)
Gene Name	MMD PAQR11
Cellular localization	Late endosome membrane ; Multi-pass membrane protein . Lysosome membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	23531
Human Swiss-Prot Number	Q15546
Alternative Names	
Background	

