

# TTL11 rabbit pAb

Cat No.:ES19911

For research use only

## Overview

Product Name	TTL11 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TTL11 AA range: 158-208
Specificity	This antibody detects endogenous levels of human TTL11
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Tubulin polyglutamylase TTLL11 (EC 6.-.-.-) (Tubulin--tyrosine ligase-like protein 11)
Gene Name	TTLL11 C9orf20
Cellular localization	Cytoplasm, cytoskeleton, cilium basal body . Cytoplasm, cytoskeleton .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	158135
Human Swiss-Prot Number	Q8NHH1
Alternative Names	
Background	

