

# PGAM1 rabbit pAb

Cat No.:ES19916

For research use only

## Overview

Product Name	PGAM1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PGAM1 AA range: 14-64
Specificity	This antibody detects endogenous levels of human PGAM1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Phosphoglycerate mutase 1 (EC 3.1.3.13) (EC 5.4.2.1) (EC 5.4.2.4) (BPG-dependent PGAM 1) (Phosphoglycerate mutase isozyme B) (PGAM-B)
Gene Name	PGAM1 PGAMA CDABP0006
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	5223
Human Swiss-Prot Number	P18669
Alternative Names	
Background	

