

# DNM3L rabbit pAb

Cat No.:ES19918

For research use only

## Overview

Product Name	DNM3L rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human DNM3L AA range: 52-102
Specificity	This antibody detects endogenous levels of human DNM3L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃. Avoid repeated freeze-thaw cycles.
Protein Name	DNA (cytosine-5)-methyltransferase 3-like
Gene Name	DNMT3L
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	29947
Human Swiss-Prot Number	Q9UJW3
Alternative Names	
Background	

