

TNPO2 rabbit pAb

Cat No.:ES19919

For research use only

Overview

Product Name	TNPO2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TNPO2 AA range: 306-356
Specificity	This antibody detects endogenous levels of human TNPO2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Transportin-2 (Karyopherin beta-2b)
Gene Name	TNPO2
Cellular localization	Cytoplasm . Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	30000
Human Swiss-Prot Number	O14787
Alternative Names	
Background	

