

C144A rabbit pAb

Cat No.:ES19930

For research use only

Overview

Product Name	C144A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human C144A AA range: 1116-1166
Specificity	This antibody detects endogenous levels of human C144A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Coiled-coil domain-containing protein 144A
Gene Name	CCDC144A KIAA0565
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	9720
Human Swiss-Prot Number	A2RUR9
Alternative Names	
Background	

