

PPC1A rabbit pAb

Cat No.:ES19943

For research use only

Overview

Product Name	PPC1A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PPC1A AA range: 165-215
Specificity	This antibody detects endogenous levels of human PPC1A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Phosphatidate phosphatase PPAPDC1A (EC 3.1.3.4) (Phosphatidic acid phosphatase type 2 domain-containing protein 1A)
Gene Name	PPAPDC1A DPPL2 PPAPDC1
Cellular localization	Membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	196051
Human Swiss-Prot Number	Q5VZY2
Alternative Names	
Background	

