

NDST3 rabbit pAb

Cat No.:ES19945

For research use only

Overview

Product Name	NDST3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NDST3 AA range: 132-182
Specificity	This antibody detects endogenous levels of human NDST3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Bifunctional heparan sulfate N-deacetylase/N-sulfotransferase 3 (EC 2.8.2.8) (Glucosaminyl N-deacetylase/N-sulfotransferase 3) (NDST-3) (hNDST-3) (N-heparan sulfate sulfotransferase 3) (N-HSST 3) [Inc
Gene Name	NDST3 HSST3 UNQ2544/PRO4998
Cellular localization	Golgi apparatus membrane ; Single-pass type II membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	9348
Human Swiss-Prot Number	O95803
Alternative Names	
Background	

