

CATF (Cleaved-Ala271) rabbit pAb

Cat No.:ES19959

For research use only

Overview

Product Name	CATF (Cleaved-Ala271) rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human CATF
	(Cleaved-Ala271)
Specificity	This antibody detects endogenous levels of
	Human,Mouse,Rat CATF (Cleaved-Ala271, protein
	was cleaved amino acid sequence between
	270-271)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20 $^\circ\!\mathrm{C}$. Avoid repeated freeze-thaw cycles.
Protein Name	CATF (Cleaved-Ala271)
Gene Name	CTSF
Cellular localization	Lysosome.
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	23 55kD
Human Gene ID	8722
Human Swiss-Prot Number	Q9UBX1
Alternative Names	Cathepsin F (CATSF;EC 3.4.22.41)
Background	



ELKbio@ELKbiotech.com

.com www

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C