



CATO (Cleaved-Leu108) rabbit pAb

Cat No.:ES19962

For research use only

Overview

Product Name	CATO (Cleaved-Leu108) rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human CATO (Cleaved-Leu108)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat CATO (Cleaved-Leu108, protein was cleaved amino acid sequence between 107-108)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	CATO (Cleaved-Leu108)
Gene Name	CTSO CTSO1
Cellular localization	Lysosome.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	23 35kD
Human Gene ID	1519
Human Swiss-Prot Number	P43234
Alternative Names	Cathepsin O (EC 3.4.22.42)
Background	

