

## Kallikrein-1 (Cleaved-Ile25) rabbit pAb

## Cat No.:ES20015

For research use only

## Overview

Product Name	Kallikrein-1 (Cleaved-Ile25) rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human;Mouse;Rat
<b>Recommended dilutions</b>	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human
	Kallikrein-1 (Cleaved-Ile25)
Specificity	This antibody detects endogenous levels of
	Human, Mouse, Rat Kallikrein-1 (Cleaved-Ile25,
	protein was cleaved amino acid sequence between
	24-25)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20 $^\circ\!\mathrm{C}$ . Avoid repeated freeze-thaw cycles.
Protein Name	Kallikrein-1 (Cleaved-Ile25)
Gene Name	KLK1
Cellular localization	
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	25kD
Human Gene ID	3816
Human Swiss-Prot Number	P06870/P20151/P07288
Alternative Names	Kallikrein-1 (EC 3.4.21.35;Kidney/pancreas/salivary
	gland kallikrein;Tissue kallikrein)
Background	Kallikreins are a subgroup of serine proteases having
	diverse physiological functions. Growing evidence
	suggests that many kallikreins are implicated in
	carcinogenesis and some have potential as novel
	cancer and other disease biomarkers. This gene is
	one of the fifteen kallikrein subfamily members



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



located in a cluster on chromosome 19. This protein is functionally conserved in its capacity to release the vasoactive peptide, Lys-bradykinin, from low molecular weight kininogen. [provided by RefSeq, Jul 2008],



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C