



# KLK1/2/3 (Cleaved-Ile25) rabbit pAb

Cat No.:ES20017

For research use only

## Overview

<b>Product Name</b>	KLK1/2/3 (Cleaved-Ile25) rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	WB; ELISA
<b>Species Cross-Reactivity</b>	Human;Rat;Mouse;
<b>Recommended dilutions</b>	WB 1:1000-2000 ELISA 1:5000-20000
<b>Immunogen</b>	Synthesized peptide derived from human KLK1/2/3 (Cleaved-Ile25)
<b>Specificity</b>	This antibody detects endogenous levels of Human KLK1/2/3 (Cleaved-Ile25, protein was cleaved amino acid sequence between 24-25 )
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at -20°C. Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	KLK1/2/3 (Cleaved-Ile25)
<b>Gene Name</b>	KLK1/2/3
<b>Cellular localization</b>	
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	25kD
<b>Human Gene ID</b>	3817
<b>Human Swiss-Prot Number</b>	P20151/P06870/P07288
<b>Alternative Names</b>	Kallikrein-2 (EC 3.4.21.35;Glandular kallikrein-1;hGK-1;Tissue kallikrein-2)
<b>Background</b>	

