

Kallikrein-6 (Cleaved-Leu22) rabbit pAb

Cat No.:ES20019

For research use only

Overview

Product Name	Kallikrein-6 (Cleaved-Leu22) rabbit pAb
Host species	Rabbit
Applications	WB; ELISA;IHC
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Immunogen	Synthesized peptide derived from human
<u> </u>	Kallikrein-6 (Cleaved-Leu22)
Specificity	This antibody detects endogenous levels of Human
	Kallikrein-6 (Cleaved-Leu22, protein was cleaved
	amino acid sequence between 21-22)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Kallikrein-6 (Cleaved-Leu22)
Gene Name	KLK6 PRSS18 PRSS9
Cellular localization	Secreted. Nucleus, nucleolus. Cytoplasm.
	Mitochondrion. Microsome. In brain, detected in
	the nucleus of glial cells and in the nucleus and
	cytoplasm of neurons. Detected in the
	mitochondrial and microsomal fractions of HEK-293
	cells and released into the
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	25kD
Human Gene ID	5653
Human Swiss-Prot Number	Q92876
Alternative Names	Kallikrein-6 (EC 3.4.21;Neurosin;Protease
	M;SP59;Serine protease 18;Serine protease 9;Zyme)
Background	

Background



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Immunohistochemical analysis of paraffin-embedded human oophoroma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C