

## MCFD2 (Cleaved-Glu27) rabbit pAb

## Cat No.:ES20034

For research use only

## Overview

Product Name	MCFD2 (Cleaved-Glu27) rabbit pAb
Host species	Rabbit
Applications	WB;IHC
Species Cross-Reactivity	Human;Mouse;Rat
<b>Recommended dilutions</b>	WB 1:500-2000;IHC-p 1:50-300
Immunogen	Synthesized peptide derived from human MCFD2
-	(Cleaved-Glu27)
Specificity	This antibody detects endogenous levels of
	Human,Mouse,Rat MCFD2 (Cleaved-Glu27, protein
	was cleaved amino acid sequence between 26-27)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	MCFD2 (Cleaved-Glu27)
Gene Name	MCFD2 SDNSF
<b>Cellular localization</b>	Endoplasmic reticulum-Golgi intermediate
	compartment . Endoplasmic reticulum . Golgi
	apparatus .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	14 16kD
Human Gene ID	90411
Human Swiss-Prot Number	Q8NI22
Alternative Names	Multiple coagulation factor deficiency protein 2
	(Neural stem cell-derived neuronal survival protein)
Background	

Background



www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C





Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C