



NOTCH3 (Cleaved-Val1662) rabbit pAb

Cat No.:ES20040

For research use only

Overview

Product Name	NOTCH3 (Cleaved-Val1662) rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human NOTCH3 (Cleaved-Val1662)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat NOTCH3 (Cleaved-Val1662, protein was cleaved amino acid sequence between 1661-1662)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	NOTCH3 (Cleaved-Val1662)
Gene Name	NOTCH3
Cellular localization	Cell membrane ; Single-pass type I membrane protein.; [Notch 3 intracellular domain]: Nucleus. Following proteolytical processing NICD is translocated to the nucleus.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	75 260kD
Human Gene ID	4854
Human Swiss-Prot Number	Q9UM47
Alternative Names	Neurogenic locus notch homolog protein 3 (Notch 3) [Cleaved into: Notch 3 extracellular truncation; Notch 3 intracellular domain]

Background

