

PRSS8 (heavy chain, Cleaved-Ile45) rabbit

pAb

Cat No.:ES20043

For research use only

Overview

Product Name	PRSS8 (heavy chain, Cleaved-Ile45) rabbit pAb
Host species	Rabbit
Applications	WB;IHC
Species Cross-Reactivity	Human; Mouse; Rat
Recommended dilutions	WB 1:500-2000;IHC-p 1:50-300
Immunogen	Synthesized peptide derived from human PRSS8
	(heavy chain, Cleaved-Ile45)
Specificity	This antibody detects endogenous levels of
	Human,Mouse,Rat PRSS8 (heavy chain,
	Cleaved-Ile45, protein was cleaved amino acid
	sequence between 44-45)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	PRSS8 (heavy chain, Cleaved-Ile45)
Gene Name	PRSS8
Cellular localization	[Prostasin]: Cell membrane; Single-pass membrane
	protein.; [Prostasin light chain]: Secreted,
	extracellular space. Found in the seminal fluid.
	Secreted after cleavage of its C-terminus.; [Prostasin
	heavy chain]: Secreted, extracellular space. Found in
	the
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	32kD
Human Gene ID	5652
Human Swiss-Prot Number	Q16651
Alternative Names	Prostasin (EC 3.4.21;Channel-activating protease
•	1;CAP1;Serine protease 8) [Cleaved into: Prostasin



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Background

light chain; Prostasin heavy chain]

Immunohistochemical analysis of paraffin-embedded human Colon cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C