



# Histone H3 (Tri Methyl Lys27) rabbit pAb

Cat No.:ES20046

For research use only

## Overview

<b>Product Name</b>	Histone H3 (Tri Methyl Lys27) rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	WB; ELISA
<b>Species Cross-Reactivity</b>	Human;Mouse;Rat
<b>Recommended dilutions</b>	WB 1:1000-2000 ELISA 1:5000-20000
<b>Immunogen</b>	Synthesized peptide derived from human Histone H3 (Tri Methyl Lys27)
<b>Specificity</b>	This antibody detects endogenous levels of Human,Mouse,Rat Histone H3 (Tri Methyl Lys27)
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at -20°C . Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	Histone H3 (Tri Methyl Lys27)
<b>Gene Name</b>	H3
<b>Cellular localization</b>	Nucleus. Chromosome.
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	17kD
<b>Human Gene ID</b>	8350
<b>Human Swiss-Prot Number</b>	P68431/Q71DI3/P84243/Q6NXT2
<b>Alternative Names</b>	Histone H3.1 (Histone H3/a;Histone H3/b;Histone H3/c;Histone H3/d;Histone H3/f;Histone H3/h;Histone H3/i;Histone H3/j;Histone H3/k;Histone H3/l)

## Background

