



Histone H3 (Di Methyl Lys37) rabbit pAb

Cat No.:ES20047

For research use only

Overview

Product Name	Histone H3 (Di Methyl Lys37) rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human Histone H3 (Di Methyl Lys37)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat Histone H3 (Di Methyl Lys37)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Histone H3 (Di Methyl Lys37)
Gene Name	H3
Cellular localization	Nucleus. Chromosome.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	17kD
Human Gene ID	8350
Human Swiss-Prot Number	P68431/Q71DI3/P84243/Q6NXT2
Alternative Names	Histone H3.1 (Histone H3/a;Histone H3/b;Histone H3/c;Histone H3/d;Histone H3/f;Histone H3/h;Histone H3/i;Histone H3/j;Histone H3/k;Histone H3/l)

Background

