

Histone H4 (Tri Methyl Lys21) rabbit pAb

Cat No.:ES20050

For research use only

Overview

(Tri Methyl Lys21) Specificity This antibody detects end Human,Mouse,Rat Histor	321) rabbit pAb
ApplicationsWB; ELISASpecies Cross-ReactivityHuman;Mouse;RatRecommended dilutionsWB 1:1000-2000 ELISA 1:ImmunogenSynthesized peptide derive (Tri Methyl Lys21)SpecificityThis antibody detects end Human,Mouse,Rat Histore	
Species Cross-ReactivityHuman;Mouse;RatRecommended dilutionsWB 1:1000-2000 ELISA 1:ImmunogenSynthesized peptide deriv (Tri Methyl Lys21)SpecificityThis antibody detects end Human,Mouse,Rat Histor	
Recommended dilutionsWB 1:1000-2000 ELISA 1:ImmunogenSynthesized peptide derive (Tri Methyl Lys21)SpecificityThis antibody detects end Human,Mouse,Rat Histor	
ImmunogenSynthesized peptide deriv (Tri Methyl Lys21)SpecificityThis antibody detects end Human,Mouse,Rat Histor	
(Tri Methyl Lys21) Specificity This antibody detects end Human,Mouse,Rat Histor	5000-20000
Specificity This antibody detects end Human,Mouse,Rat Histor	ed from human Histone H4
Human, Mouse, Rat Histor	
	ogenous levels of
	e H4 (Tri Methyl Lys21)
Formulation Liquid in PBS containing 5	0% glycerol, 0.5% BSA and
0.02% sodium azide.	
Storage Store at -20°C. Avoid repe	eated freeze-thaw cycles.
Protein Name Histone H4 (Tri Methyl Lys	\$21)
Gene Name HIST1H4A H4/A H4FA; HIS	5T1H4B H4/I H4FI;
HIST1H4C H4/G H4FG; HIS	ST1H4D H4/B H4FB;
HIST1H4E H4/J H4FJ; HIST	1H4F H4/C H4FC;
HIST1H4H H4/H H4FH; HI	ST1H4I H4/M H4FM;
HIST1H4J H4/E H4FE; HIS	「1H4K H4/D H4FD;
HIST1H4L H4/K H4FK; HIS	T2H4A H4/N H4F2 H4FN
HIST2H4; HI	
Cellular localization Nucleus. Chromosome.	
Purification The antibody was affinity	purified from rabbit
antiserum by affinity-chro	matography using
epitope-specific immuno	gen.
Clonality Polyclonal	
Concentration 1 mg/ml	
Observed band 12kD	
Human Gene ID 121504	
Human Swiss-Prot Number P62805	
Alternative Names Histone H4	
Background function:Core component	of nucleosome.
Nucleosomes wrap and co	omnact DNA into
chromatin, limiting DNA a	



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling., PTM: Acetylation at Lys-6, Lys-9, Lys-13 and Lys-17 occurs in coding regions of the genome but not in heterochromatin., PTM: Citrullination at Arg-4 by PADI4 impairs methylation., PTM: Monomethylated, dimethylated or trimethylated at Lys-21. Monomethylation is performed by SET8. Trimethylation is performed by SUV420H1 and SUV420H2 and induces gene silencing., PTM: Monomethylation at Arg-4 by PRMT1 favors acetylation at Lys-9 and Lys-13. Demethylation is performed by JMJD6., PTM: Sumoylated, which is associated with transcriptional repression., PTM: Ubiquitinated by the CUL4-DDB-RBX1 complex in response to ultraviolet irradiation. This may weaken the interaction between histones and DNA and facilitate DNA accessibility to repair proteins., similarity: Belongs to the histone H4 family., subunit: The nucleosome is a histone octamer containing two molecules each of H2A, H2B, H3 and H4 assembled in one H3-H4 heterotetramer and two H2A-H2B heterodimers. The octamer wraps approximately 147 bp of DNA.,



ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C