

ERα (Acetyl Lys268) rabbit pAb

Cat No.:ES20064

For research use only

Overview

Product Name	ERα (Acetyl Lys268) rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human ER $lpha$
<u> </u>	(Acetyl Lys268)
Specificity	This antibody detects endogenous levels of Human
	ERa (Acetyl Lys268)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20 $^\circ\!\mathrm{C}$. Avoid repeated freeze-thaw cycles.
Protein Name	ERα (Acetyl Lys268)
Gene Name	ESR1 ESR NR3A1
Cellular localization	[Isoform 1]: Nucleus . Cytoplasm . Cell membrane ;
	Peripheral membrane protein ; Cytoplasmic side . A
	minor fraction is associated with the inner
	membrane.; [Isoform 3]: Nucleus. Cytoplasm. Cell
	membrane; Peripheral membrane protein;
	Cytoplasmic side. Cell membrane; Single-pass type I
	membrane protein. Associated with the inner
	membrane via palmitoylation (Probable). At least a
	subset exists as a transmembrane protein with a
	N-terminal extracellular domain; Nucleus. Golgi
	apparatus. Cell membrane. Colocalizes with ZDHHC7
	and ZDHHC21 in the Golgi apparatus where most
	probably palmitoylation occurs. Associated with the
	plasma membrane when palmitoylated.
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



Observed band 65kD Human Gene ID 2099 Human Swiss-Prot Number P03372 Estrogen receptor (ER;ER-alpha;Estradiol Alternative Names receptor;Nuclear receptor subfamily 3 group A member 1) Background domain:Composed of three domains: a modulating N-terminal domain, a DNA-binding domain and a C-terminal steroid-binding domain., function: Nuclear hormone receptor. The steroid hormones and their receptors are involved in the regulation of eukaryotic gene expression and affect cellular proliferation and differentiation in target tissues., online information: Estrogen receptor entry, polymorphism: Genetic variations in ESR1 are correlated with bone mineral density (BMD). Low BMD is a risk factor for osteoporotic fracture. Osteoporosis is characterized by reduced bone mineral density, disrutption of bone microarchitecture, and the alteration of the amount and variety of non-collagenous proteins in bone. Osteoporotic bones are more at risk of fracture., PTM: Glycosylated; contains N-acetylglucosamine, probably O-linked., PTM: Phosphorylated by cyclin A/CDK2. Phosphorylation probably enhances transcriptional activity., similarity: Belongs to the nuclear hormone receptor family., similarity: Belongs to the nuclear hormone receptor family. NR3 subfamily., similarity: Contains 1 nuclear receptor DNA-binding domain., subunit: Interacts with SLC30A9 (By similarity). Binds DNA as a homodimer. Can form a heterodimer with ESR2. Interacts with NCOA3, NCOA5 and NCOA6 coactivators, leading to a strong increase of transcription of target genes. Interacts with NCOA7 in a ligand-inducible manner. Interacts with PHB2, PELP1 and UBE1C. Interacts with AKAP13. Interacts with CUEDC2. Interacts with KDM5A. Interacts with SMARD1. Interacts with HEXIM1 and MAP1S. Interacts with PBXIP1. Interaction with MUC1 is stimulated by 7 beta-estradiol (E2) and enhances ERS1-mediated



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.G



transcription. Interacts with DNTTIP2, FAM120B and UIMC1. Interacts with isoform 4 of TXNRD1. Interacts with MLL2. Interacts with ATAD2 and this interaction is enhanced by estradiol.,



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C