

Histone H2A.Z (Acetyl Lys12) rabbit pAb

Cat No.:ES20073

For research use only

Overview

Product Name	Histone H2A.Z (Acetyl Lys12) rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human Histone
	H2A.Z (Acetyl Lys12)
Specificity	This antibody detects endogenous levels of
	Human,Mouse,Rat Histone H2A.Z (Acetyl Lys12)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20 $^\circ\!\mathrm{C}$. Avoid repeated freeze-thaw cycles.
Protein Name	Histone H2A.Z (Acetyl Lys12)
Gene Name	H2AFZ H2AZ
Cellular localization	Nucleus. Chromosome.
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	15kD
Human Gene ID	3015
Human Swiss-Prot Number	P0C0S5
Alternative Names	Histone H2A.Z (H2A/z)
Background	function:Variant histone H2A which replaces
	conventional H2A in a subset of nucleosomes.
	Nucleosomes wrap and compact DNA into
	chromatin, limiting DNA accessibility to the cellular
	machineries which require DNA as a template.
	Histones thereby play a central role in transcription
	regulation, DNA repair, DNA replication and
	chromosomal stability. DNA accessibility is regulated
	via a complex set of post-translational modifications



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



of histones, also called histone code, and nucleosome remodeling. May be involved in the formation of constitutive heterochromatin. May be required for chromosome segregation during cell division., mass spectrometry: Monoisotopic, not modified PubMed:16457589,PTM:Acetylated on Lys-5, Lys-8 and Lys-12 during interphase. Acetylation disappears at mitosis., PTM: Monoubiquitination of Lys-122 gives a specific tag for epigenetic transcriptional repression.,PTM:Not phosphorylated., similarity: Belongs to the histone H2A family., subunit: The nucleosome is a histone octamer containing two molecules each of H2A, H2B, H3 and H4 assembled in one H3-H4 heterotetramer and two H2A-H2B heterodimers. The octamer wraps approximately 147 bp of DNA. H2A or its variant H2AFZ forms an heterodimer with H2B. H2AFZ interacts with INCENP.,



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C