



IL27B rabbit pAb

Cat No.:ES20343

For research use only

Overview

Product Name	IL27B rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Mouse
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from mouse IL27B
Specificity	This antibody detects endogenous levels of Mouse IL27B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	IL27B
Gene Name	IL27B
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	O35228
Alternative Names	Interleukin-27 subunit beta (IL-27 subunit beta;IL-27B;Epstein-Barr virus-induced gene 3 protein;EBV-induced gene 3 protein)

Background

