



CXL10 rabbit pAb

Cat No.:ES20365

For research use only

Overview

Product Name	CXL10 rabbit pAb
Host species	Rabbit
Applications	WB; ELISA
Species Cross-Reactivity	Mouse
Recommended dilutions	WB 1:1000-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from mouse CXL10
Specificity	This antibody detects endogenous levels of Mouse CXL10
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	CXL10
Gene Name	CXL10
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	
Alternative Names	C-X-C motif chemokine 10 (10 kDa interferon gamma-induced protein;Gamma-IP10;IP-10;Small-inducible cytokine B10) [Cleaved into: CXCL10(1-73)]

Background

