



Histone H2A (Acetyl Lys15) rabbit pAb

Cat No.:ES20537

For research use only

Overview

Product Name	Histone H2A (Acetyl Lys15) rabbit pAb
Host species	Rabbit
Applications	WB
Species	Human;Mouse;Rat
Cross-Reactivity	
Recommended dilutions	WB: 1:500-1000
Immunogen	Synthetic Peptide of Histone H2A (Acetyl Lys15)
Specificity	The antibody detects endogenous Histone H2A (Acetyl Lys15) protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Histone H2A type 1/Histone H2A type 2/Histone H2A type 3
Gene Name	HIST1H2AG/HIST1H2AI/HIST1H2AK/HIST1H2AL/HIST1H2AM/HIST2H2AA3/HIST2H2AA4/HIST3H2A
Cellular localization	Nucleus. Chromosome.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	
Observed band	14kD
Human Gene ID	8329/8330/8332/8336/8969/723790/8337/92815





Human POC0S8/Q6FI13/Q7L7L0

**Swiss-Prot
Number**

**Alternative
Names** HIST1H2AG; H2AFP; HIST1H2AI; H2AFC; HIST1H2AK; H2AFD;
HIST1H2AL; H2AFI; HIST1H2AM; H2AFN; Histone H2A type 1; H2A.1;
Histone H2A/p; HIST2H2AA3; H2AFO; HIST2H2AA; HIST2H2AA4;
Histone H2A type 2-A; Histone H2A.2; Histone H2A/o; HIST3H2A;
Histone H2A type 3;

Background Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H2A family. Transcripts from this gene lack polyA tails but instead contain a palindromic termination element. This gene is found in the small histone gene cluster on chromosome 6p22-p21.3. [provided by RefSeq, Aug 2015],

Western blot analysis of extracts from HeLa cells, untreated (-) or treated, 1:5000. Secondary antibody(catalog#:RS0002) was diluted at 1:20000

