

GAL4 DNA Binding Domain rabbit pAb

Cat No.:ES20743

For research use only

Overview

Product Name	GAL4 DNA Binding Domain rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	
Recommended dilutions	WB: 1:10000-50000
Immunogen	Recombinant Protein of GAL4 DNA Binding Domain
Specificity	The antibody detects endogenous GAL4 DNA
	Binding Domain protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	
Gene Name	
Cellular localization	
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	
Observed band	
Human Gene ID	
Human Swiss-Prot Number	
Alternative Names	YM3411
Background	The GAL4 protein of Saccharomyces cerevisiae is
	one of the most thoroughly characterized
	transcriptional activators. Since the N-terminal 147
	amino acid residues of GAL4 are sufficient to
	mediate specific and strong binding to DNA, but are
	incapable of efficient transcriptional activation , this
	protein fragment has frequently been used to confer
	specific DNA binding in experiments examining
	transcriptional activation functions of heterologous
	proteins. This approach is facilitated by the finding
3	



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



that higher eukaryotes lack endogenous proteins that enhance transcription from the consensus GAL4-binding site. Fusions between GAL4 (aa 1-147) and activating domains from a variety of transcriptional regulatory proteins can activate transcription in yeast, plant, insects and mammalian cells. A unique "two-hybrid" system has been developed using GAL4 fusions in yeast to identify specific protein-protein interactions.

Western blot analysis of Recombinant BD protein using GAL4 DNA Binding Domain Polyclonal Antibody. Secondary antibody(catalog#:RS0002) was diluted at 1:20000





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C