

CYB561D2 rabbit pAb

Cat No.: ES4588

For research use only

Overview

Purification

Product Name CYB561D2 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human;Rat;Mouse;

Recommended dilutions Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not

yet tested in other applications.

Immunogen The antiserum was produced against synthesized

peptide derived from human C56D2. AA

range:55-104

Specificity CYB561D2 Polyclonal Antibody detects endogenous

levels of CYB561D2 protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Cytochrome b561 domain-containing protein 2

Gene Name CYB561D2

Cellular localization Endoplasmic reticulum membrane ; Multi-pass

membrane protein . Cytoplasmic vesicle membrane ; Multi-pass membrane protein . The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 24kD
Human Gene ID 11068
Human Swiss-Prot Number O14569

Alternative Names CYB561D2; 101F6; LUCA12.2; Cytochrome b561

domain-containing protein 2; Putative tumor

suppressor protein 101F6

Background cofactor:Binds 2 heme groups


non-covalently., similarity: Contains 1 cytochrome

b561 domain.,




+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com





Western Blot analysis of Jurkat cells using CYB561D2 Polyclonal Antibody



Western blot analysis of lysates from Jurkat cells, using C56D2 Antibody. The lane on the right is blocked with the synthesized peptide.

