

IL-27A rabbit pAb

Cat No.:ES8732

For research use only

Overview

Product Name	IL-27A rabbit pAb
Host species	Rabbit
Applications	IHC;IF;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	IHC-p 1:50-200, ELISA 1:10000-20000
Immunogen	Synthetic peptide from human protein at AA range:
	1-50
Specificity	The antibody detects endogenous IL-27A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Interleukin-27 subunit alpha (IL-27 subunit alpha)
	(IL-27-A) (IL27-A) (p28)
Gene Name	IL27 IL27A
Cellular localization	Secreted . Does not seem to be secreted without
	coexpression of EBI3.
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	246778
Human Swiss-Prot Number	Q8NEV9
Alternative Names	Interleukin-27 subunit alpha (IL-27 subunit
	alpha;IL-27-A;IL27-A;p28)
Background	The protein encoded by this gene is one of the
	subunits of a heterodimeric cytokine complex. This
	protein is related to interleukin 12A (IL12A). It
	interacts with Epstein-Barr virus induced gene 3
	(EBI3), a protein similar to interleukin 12B (IL12B),
	and forms a complex that has been shown to drive
	rapid expansion of naive but not memory CD4(+) T



+86-27-59760950

ELKbio@ELKbiotech.com


www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



cells. The complex is also found to synergize strongly with interleukin 12 to trigger interferon gamma (IFNG) production of naive CD4(+) T cells. The biological effects of this cytokine are mediated by the class I cytokine receptor (WSX1/TCRR). [provided by RefSeq, Jul 2008],

Immunohistochemical analysis of paraffin-embedded human-lung, antibody was diluted at 1:200





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C