



Tmem72 rabbit pAb

Cat No.:ES8945

For research use only

Overview


Product Name	Tmem72 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:500-2000, ELISA 1:10000-20000
Immunogen	Synthesized peptide derived from human Tmem72 Polyclonal
Specificity	This antibody detects endogenous levels of Tmem72.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Transmembrane protein 72 (Kidney-specific secretory protein of 37 kDa)
Gene Name	TMEM72 C10orf127 KSP37
Cellular localization	Membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	30kD
Human Gene ID	643236
Human Swiss-Prot Number	A0PK05
Alternative Names	Transmembrane protein 72 (Kidney-specific secretory protein of 37 kDa)
Background	This gene encodes a transmembrane protein which may be expressed specifically in the kidney. [provided by RefSeq, Sep 2016],





ELK Biotechnology

Western blot analysis of various lysate, antibody was diluted at 1000. Secondary antibody(catalog#:RS0002) was diluted at 1:20000



1 mouse-liver

2 CACO2



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C