



LRP12 rabbit pAb

Cat No.:ES9776

For research use only

Overview

Product Name	LRP12 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 90-170
Specificity	LRP12 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Low-density lipoprotein receptor-related protein 12 (LRP-12) (Suppressor of tumorigenicity 7 protein)
Gene Name	LRP12 ST7
Cellular localization	Membrane ; Single-pass type I membrane protein . Membrane, coated pit .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	94kD
Human Gene ID	29967
Human Swiss-Prot Number	Q9Y561
Alternative Names	
Background	This gene encodes a member of the low-density lipoprotein receptor related protein family. The product of this gene is a transmembrane protein that is differentially expressed in many cancer cells. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Feb 2010],

