

## **Recombinant Human VEGF165**

Catalog #	EPT029
Expression Host	Human Cells
DESCRIPTION	Recombinant Human Vascular Endothelial Growth
	Factor A is produced by our Mammalian expression
	system and the target gene encoding Ala27-Arg191 is
	expressed.
Accession	P15692-4
Synonyms	Vascular Endothelial Growth Factor Isoform 165;
	VEGF165
Mol Mass	19.1 KDa
AP Mol Mass	18-22 KDa, reducing conditions
Purity	Greater than 95% as determined by reducing
	SDS-PAGE.
Endotoxin	Less than 0.001 ng/ $\mu$ g (0.01 EU/ $\mu$ g) as determined by
	LAL test.
FORMULATION	Lyophilized from a 0.2 $\mu\text{m}$ filtered solution of 20mM
	PB, 150mM NaCl, pH 7.2.
RECONSTITUTION	Always centrifuge tubes before opening.Do not mix by



+86-27-59760950 ELKbio@ELKbiotech.com www.elkbiotech.com 23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



vortex or pipetting.

It is not recommended to reconstitute to a concentration less than 100µg/ml.

Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SHIPPING The product is shipped at ambient temperature.Upon receipt, store it immediately at the temperature listed below.

STORAGELyophilized protein should be stored at < -20 ° C,<br/>though stable at room temperature for 3 weeks.<br/>Reconstituted protein solution can be stored at 4-7°C<br/>for 2-7 days.

Aliquots of reconstituted samples are stable at < -20° C for 3 months.

BACKGROUND Human Vascular endothelial growth factor (VEGF), also known as VEGF-A and vascular permeability factor (VPF), belongs to the platelet-derived growth factor family of cysteine-knot growth factors. It is a potent activator in vasculogenesis and angiogenesis both physiologically and pathologically. VEGF-A has 8 differently spliced isoforms, of which VEGF165 is the




+86-27-59760950 ELKbio@ELKbiotech.com

om www.e

www.elkbiotech.com



most abundant one. VEGF165 is a disulfide-linked homodimer consisting of two glycosylated 165 amino acid polypeptide chains. VEGF stimulates the cellular response through binding to tyrosine kinase receptors VEGFR1 and VEGFR2 on the cell surface. It is widely accepted that VEGFR2 mediate almost all of the known cellular responses to VEGF while the function of VEGFR1 is less defined and is thought to modulate the VEGFR2 signaling.



SDS-PAGE



ELKbio@ELKbiotech.com +86-27-59760950

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.Q