

Recombinant Human TNF alpha

Catalog # EPT038

Expression Host E.coli

DESCRIPTION Recombinant Human Tumor Necrosis Factor Alpha is

produced by our E.coli expression system and the

target gene encoding Val77-Leu233 is expressed.

Accession P01375

Synonyms Tumor Necrosis Factor; Cachectin; TNF-Alpha; Tumor

Necrosis Factor Ligand Superfamily Member 2; TNF-a;

TNF; TNFA; TNFSF2

Mol Mass 17.5 KDa

AP Mol Mass 16 KDa, reducing conditions

Purity Greater than 95% as determined by reducing

SDS-PAGE.

Endotoxin Less than 0.1 ng/μg (1 EU/μg) as determined by LAL

test.

FORMULATION Lyophilized from a 0.2 µm filtered solution of 20mM

PB, 150mM NaCl, pH 7.0.

RECONSTITUTION Always centrifuge tubes before opening.Do not mix by



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com



vortex or pipetting.

It is not recommended to reconstitute to a concentration less than 100µg/ml.

Dissolve the lyophilized protein in distilled water.

Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SHIPPING

The product is shipped at ambient temperature.

Upon receipt, store it immediately at the temperature listed below.

STORAGE

Lyophilized protein should be stored at < -20 ° C, though stable at room temperature for 3 weeks.

Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

Aliquots of reconstituted samples are stable at < -20° C for 3 months.

BACKGROUND


TNF α is a homotrimer with a subunit molecular mass of 17 kD cytokine that binds to TNFRSF1A/TNFR1 and TNFRSF1B/TNFBR. It is mainly secreted by macrophages and can induce cell death of certain tumor cell lines. It plays a major role in growth regulation, differentiation, inflammation, viral replication, tumorigenesis, autoimmune diseases and





in viral, bacterial, fungal, and parasitic infections.

Besides inducing hemorrhagic necrosis of tumors, TNF was found to be involved in tumorigenesis, tumor metastasis, viral replication, septic shock, fever, inflammation, and autoimmune diseases including Crohn's disease, and rheumatoid arthritis as well as graft-versus-host disease.



SDS-PAGE

