

Recombinant Mouse CCL2

Catalog # EPT150

Expression Host E.coli

DESCRIPTION Recombinant Mouse C-C Motif Chemokine 2 is

produced by our E.coli expression system and the

target gene encoding Gln24-Arg96 is expressed.

Accession P10148

Synonyms C-C motif chemokine 2; Monocyte chemoattractant

protein 1; Monocyte chemotactic protein 1;

MCP-1;Platelet-derived growth factor-inducible

protein JE; Small-inducible cytokine A2; Ccl2; Je; Mcp1;

Scya2

Mol Mass 8.5 KDa

AP Mol Mass 12 KDa, reducing conditions

Purity Greater than 95% as determined by reducing

SDS-PAGE.

Endotoxin Less than 0.001 ng/ μ g (0.01 EU/ μ g) as determined by

LAL test.

FORMULATION Lyophilized from a 0.2 µm filtered solution of PBS, pH



www.elkbiotech.com



7.4.

RECONSTITUTION

Always centrifuge tubes before opening.Do not mix by vortex or pipetting.

It is not recommended to reconstitute to a concentration less than 100µg/ml.

Dissolve the lyophilized protein in distilled water.

Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SHIPPING

The product is shipped at ambient temperature.

Upon receipt, store it immediately at the temperature listed below.

STORAGE

Lyophilized protein should be stored at < -20 ° C, though stable at room temperature for 3 weeks.

Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

Aliquots of reconstituted samples are stable at < -20° C for 3 months.


BACKGROUND

C-C motif chemokine 2 (CCL2) is a member of the C-C or β chemokine family. Mouse CCL2 shares 82% amino acid (aa) identity with rat CCL2 over the entire sequence, and 58%, 56%, 55%, 53% and 53% aa identity with human, equine, porcine, bovine and





canine CCL2, respectively. Fibroblasts, glioma cells, smooth muscle cells, endothelial cells, lymphocytes and mononuclear phagocytes can produce CCL2 either constitutively or upon mitogenic stimulation, but monocytes and macrophages appear to be the major source. In addition to its chemotactic activity, CCL2 induces enzyme and cytokine release by monocytes, NK cells and lymphocytes, and histamine release by basophils that express its receptor, CCR2. Additionally, it promotes Th2 polarization in CD4+ T cells. CCL2-mediated recruitment of monocytes to sites of inflammation is proposed to play a role in the pathology of atherosclerosis, multiple sclerosis and allergic asthma.



SDS-PAGE

+86-27-59760950



ELKbio@ELKbiotech.com